

## ROKAnol<sup>®</sup> PW

<b>CHEMICAL NAME</b>	Mixture of chemicals
<b>INCI NAME</b>	-
<b>CAS NUMBER</b>	-
<b>FUNCTION</b>	Virucidal preparation
<b>TECHNICAL REQUIREMENTS</b>	Appearance at (20±25)°C..... colorless liquid pH of product..... 6 ÷ 9
<b>GENERAL DATA</b>	Active substance ..... ethanol 75g/100g Solubility.....easily soluble in water Density at 20°C, g/mL ..... approx. 0.9
<b>APPLICATION</b>	A virucidal fluid, designed for effective disinfection of objects and surfaces resistant to alcohol. Especially recommended for use in public facilities, including hospitals, offices and other public units. The product is intended for individual use and for professional use. The product can be used on all surfaces resistant to alcohol. Alcohol sensitive substrates should be pre-tested in a small area. Do not use the product to clean acrylic glass (plexiglass), painted surfaces and stained surfaces.